



CIC & MG Corporation Partnership Project

Malaysia University, Lab Tour Photos and Coffee Leaf Sample Preparation

May 7th, 2018



Our Partnership values



To Succeed we must always be striving to know and do the right things at all times. As partners it is important that we have a set of values which encourage us to behave ethically, with integrity and mutual respect.

› Teamwork - *We 're better together*

Success comes from collaboration: getting the best out of each other and supporting one another.

› Trust - *Confidence in each other*

We earn and give trust through our actions. We make promises and keep them, are open and honest about what we do.

› Innovation - *Better every day*

We innovate every day - in our product technology, our operations, and the way we work. We create solutions for our society and customers through small changes as well as major breakthroughs.

› Precision - *Excellence always*

We relentlessly pursue excellence, are thorough and diligent in what we do and always expect the best of each other.






Sample Preparation_1



Leave Sample
(After dusting & wiping)



Leave Sample- After grind



Acid Digest - 0 min
(Added HNO_3 & H_2SO_4)




1. The leaves must be cleaned properly prior to testing
 - Remove all unwanted particles from the leaves
2. Dry the samples using oven at 103°C for minimum of 8hrs
 - Oven dry the samples before grinding.
 - The oven needs to be able to attain and hold temperature of 103°C . This is rather important for an accurate reading.
 - The time to dry is typically about 8hrs for 100g samples.
3. Grind the samples as shown above
4. Acid digestion of the leaves
 - Nitric acid is rarely used alone. It is best used in combination with sulfuric acid for sample digestion.



Sample Preparation_2



Acid Digest- 30min



Acid Digest- 60 min



Acid Digest - 120 min



Filtrate- Clear Liquid



Filtrating Digested Sample



Acid Digest- 120 min (close-up)



- Care should be exercised and the literature consulted before attempting to use nitric acid in combination with other acids or organic sample digestions.



Sample Preparation_3





Residue of Digest



Disposed



Clear Solution-
Ready for ICP
Analysis



Leaf nutrient Target Elements

- Nitrogen , Phosphorus , Potassium , Calcium, Iron,
- Boron , Zinc , Manganese , Aluminum

Lab Tour and Coffee Leaf Samples



UTM Professor Shree



Coffee Leaf Samples



Chemist checking samples



Lab Tour and Analysis Team



Analysis Equipment

